



# Mini-100 Orbital-Genie™

Compact slow orbital shaking on a non-slip platform

Non-slip, platform designed with spill-proof edges is ideal for gentle shaking applications such as Western Blots, with no additional accessories necessary. Variable shaking speed control up to 100RPM.



	100 RPM with Ultra High Torque
	Super Compact Design
	Variable Speed
	Gentle Orbital Shaking

## Accessories

A non-slip tray that can be stacked onto the Mini-100 Orbital-Genie's existing platform tray, increasing the capacity for mixing by stacking up two trays. No tools required, stays in place, simply lift-to-remove.



Add stacking trays

Double or triple the capacity of your Mini-100



## Specifications

	<b>Mini-100 Orbital-Genie</b>
<b>Speed (RPM)</b>	20 – 100 RPM
<b>Orbit</b>	19mm
<b>Weight</b>	3.5 Kg (7.8 lbs.)
<b>Dimensions</b>	350 x 350mm (14 x 14in)
<b>Platform</b>	290 x 290mm (11.4 x 11.4in)
<b>Platform Weight Capacity</b>	3 Kg (6.6 lbs)
<b>Markings</b>	CE

## Mini Orbital-Genie Ordering Information

	<b>Catalog No.</b>
<b>Voltage &amp; Plug configuration</b>	<b>Mini-100 Orbital-Genie</b>
<b>120V - U.S Standard</b>	SI-M100
<b>230V - No Plug</b>	SI-M101
<b>230V - European Plug</b>	SI-M102
<b>230V - British Plug</b>	SI-M103
<b>230V - Swiss Plug</b>	SI-M104
<b>240V - Aust. Plug</b>	SI-M106
<b>100V</b>	SI-M105

## Accessory Ordering Information

<b>Catalog No.</b>	<b>Description</b>
<b>SI-M150</b>	Stacking Tray, for the Mini-100 Orbital-Genie



80 Orville Drive,  
Suite 102,  
Bohemia,  
New York 11716 USA



Phone: (631) 567-4700  
Toll Free: (888) 850-6208



Customer Service Email:  
[custsvc@scientificindustries.com](mailto:custsvc@scientificindustries.com)



Fax: (631) 567-5896

