



# Large Sky Line Analog Orbital Shaker S-3.02L



## Product Description

The Sky Line series is a high-performance line of orbital and rocking shakers capable of non-stop operation in a wide variety of medical, chemical, and biological laboratories. Designed for flexibility, these units come with a great selection of optional accessories and attachments, and in an assortment of sizes and loading capacities to satisfy your specific lab needs. With an excellent record of reliability and long life, the Sky Line series will be a superior partner in your laboratory.

## Features

- Optimal amplitude of platform rotation for intensive mixing
- User-friendly controls
- High-precision, rotational speed controls

## Specifications

Size of Working Space of Platform	415 x 295
Amplitude of Rotation of the Platform (Rocking Angle)	10 / 20
Type of Rotation	Orbital
Speed of Rotation (Rocking)	from 1 to 500 / 300
Environmental Temperature	from +4 to +45
Relative Moisture of Atmosphere (at the temperature of 20 °C)	up to 80
Power Supply (Adapter)	DC 24V, 1A
Size of the Device (length x width x height)	415 x 365 x 100
Weight	7.7/8.7



SHKR-15

Standard Microplate Adapter



SHKR-17

Universal Roller Clip



SHKR-19

Rubber Anti-Slip Mat



SHKR-14

30mm Test Tube Adapter



SHKR-13

16mm Test Tube Adapter



SHKR-12

Magnetic Clip Box

