

## MS7-S 7" Square Magnetic Stirrer

### Features

- Speed range of 0-1500rpm
- Max. stirring quantity of H<sub>2</sub>O at 10L
- Ceramic work plate provides excellent chemical resistant performance



Chemical resistance



Strong magnetism

### Specifications

Work plate dimension[WxD]	
Work plate material	
Motor type	
Motor input power	
Motor output power	
Power	
Voltage	
Stirring positions	
Max. stirring quantity[H <sub>2</sub> O]	
Max. magnetic bar[length]	
Speed range	
Speed display	
Protection class	
Dimension[WxDxH]	
Weight	
Permissible ambient temperature and humidity	

## MS-S Magnetic Stirrer

### Features

- Maintenance free brushless DC motor
- Speed range of 0-1500 rpm
- Max. stirring quantity of H<sub>2</sub>O at 20L
- Stainless steel work plate cover with ceramic material provides good chemical resistance performance



### F100 protective cover

for non-hotplate model,  
use with MS-S



MS7-S	MS-S
184x184mm(7inch)	φ135mm(5 inch)
Glass ceramic	Stainless steel cover with ceramic
Shaded pole motor	Brushless DC motor
15W	18W
1.5W	10W
30W	30W
100-120/200-240V,50/60Hz	100-240V,50/60Hz
1	1
10L	20L
80mm	80mm
0-1500rpm	0-1500rpm
Scale	Scale
IP21	IP42
215x360x112mm	160x280x85mm
3.8kg	2.8kg
5-40°C, 80%RH	5-40°C, 80%RH