

4 Cu. Ft. Standard Undercounter Freezer Freestanding, ABT-HS-UCFS-0420W



- Standard Undercounter Freezer Freestanding
- 4 Cu. Ft. Capacity
- One solid swing door, right hinged, not reversible
- Mechanical thermostat
- Adjustable operating temperature range: -15°C – -25°C
- Four drawers
- Keyed door lock
- Magnetic door gasket for positive seal
- Recessed door handle
- Leveling legs
- Probe access port (3/8")
- White cabinet and door
- High density urethane foam cabinet and door insulation
- HFC-free refrigerant (R600a), foam insulation and packaging
- One year parts and labor warranty
- Exterior dimensions: 20"W x 22 5/8"D x 33 3/4"H
- Note: This unit must have 4" clearance on sides and back for adequate ventilation
- Manual defrost
- 115V, 60 Hz, .90 Amps, 1/7 HP
- UL/C-UL listed