

Catalog Number: 430172

Product Description: Corning® 16 x 125 mm tissue culture tube, plug seal cap

Product Dimensions:

Length w/cap	-	4.90 in.	O.D. (cap)	-	0.84 in.
Length w/o cap	-	4.86in.	Tolerences	-	+/-0.05 in.
Width (tube)	-	0.58 in.			

Sterilization:

This lot has been irradiated and dosimetrically released based on ANSI/AAMI/ISO 11137 *Sterilization of healthcare products-Requirements for validation and routine control-Radiation sterilization*.

Sterility Assurance Level: SAL 10⁻⁵

Pyrogens:

The product has been tested and has met the criteria established in the current version of ANSI/AAMI ST 72: *Bacterial Endotoxins - Test methodologies, routine testing, and alternative to batch testing*. Results: less than 0.1EU/mL.

Lot Number Designation:

8 Digit Lot Number: First 3 digits - Julian date, start of manufacturing; Next 2 digits - Year of manufacture; Last 3 digits - Batch identification.

Revision Date:

10-31-18

Rev No: 2