

-20°C MANUAL DEFROST FREEZER

Packed with enough cooling power to maintain -20° C throughout with $\pm 2^{\circ}\text{C}$ temperature uniformity, this 25 ft³ freezer has six epoxy coated shelves to optimize storage capacity and a digital microprocessor to display the temperature. The chamber is sized to slide through door frames and plugs into a standard wall outlet. A floor drain with tubing attached is included for defrosting. Options such as interior circulating fan for $\pm 0.5^{\circ}\text{C}$ uniformity, extra shelves, door lock, hi/low alarm with battery back-up or relay, temperature chart recorder, or stainless steel exterior construction are available.

- Six 28" w x 24" d shelves, epoxy coated
- Microprocessor temp control
- Digital temperature display
- -24° thru 0°C temp range, $\pm 2^{\circ}\text{C}$
- Constant temperature cold wall refrigeration system
- Casters, locking, 4"
- Powder coated white or gray steel exterior
- White coated aluminum interior
- Non-CFC foamed insulation
- Self-closing, hinged, solid foamed door
- One year warranty on parts and labor, four add'l years on compressor motor (domestic)



Model	Capacity (ft ³ /liters)	Dimensions			# Shelves and Size	Electrical (domestic)	HP	Ship Weight (lbs/kg)
		Exterior (in/cm)		Interior (in/cm)				
FMS33SD	25/708	33/84	W	29/74	6 28" x 24"	115/60/1 15 amps	1/2	400/180
		29/74	D	26/66				
		74.5/189	H	57.5/146				

Refrigerant is 404A.