Product Name:	PH-ABT-HC-23G
Product Name: Description:	Premier Pharmacy/Vaccine Glass Door Refrigerator 23 Cu. Ft. capacity One swing glass door, self-closing (energy efficient, double pane vacuum insulated), right hinged, not reversible Microprocessor Temperature Controller Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C) Temperature monitor device (TMD) complies with the current CDC guidelines, with 3 years certification of calibration, "buffered" probe in the product simulated solution, min/max memory. F/C switchable, field installable, and visual & audible temp alarm Standard Digital Temperature Display and Alarm Module with battery back-up to support temperature monitoring and alarms Sample simulated digital temperature display; .1° Resolution Adjustable temperature display in Celsius or Fahrenheit Audible and visual high/low temperature alarms 2 probes; 1 in air and 1 in sample bottle Sensor error alarm Power failure alarm Resettable 1-point Min/Max History Remote alarm contacts Alarm validation (audible, visual, and remote alarm) Alarm mute with ring back 6 total shelves, 4 adjustable plus two additional Vaccine storage power cord warning label Forced draft air circulation Keyed door lock Magnetic door gasket for positive seal Heavy duty door handle LED interior lights are safety shielded and switch controlled Swivel Casters (3 ½") Probe access port (¾") Powder coated steel white exterior and interior
	Magnetic door gasket for positive seal Heavy duty door handle LED interior lights are safety shielded and switch controlled Swivel Casters (3 <sup>1</sup> / <sub>2</sub> ") Probe access port ( <sup>3</sup> / <sub>8</sub> ")
	Hydrocarbon, natural refrigerant (R290) Pharmacy refrigerator/freezer toolkit and temperature logs Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 26 <sup>7</sup> / <sub>8</sub> " W x 34 <sup>3</sup> / <sub>4</sub> " D x 81 <sup>3</sup> / <sub>4</sub> " H with casters; 78 <sup>1</sup> / <sub>4</sub> " H without casters Cycle defrost
	115V, 60 Hz, 3 Amps, 1/5 HP ETL listed Shipping weight: 336 lbs.