Product Name:	PH-ABT-HC-12G
Description:	Premier Pharmacy/Vaccine Glass Door Refrigerator
	12 Cu. Ft. capacity
	One swing glass door, self-closing (energy efficient, double pane
	vacuum insulated), right hinged, not reversible
	Microprocessor Temperature Controller
	Adjustable operating temperature range: 36°F – 46°F (2°C – 8°C)
	Temperature monitor device (TMD) complies with the current CDC
	guidelines, with 3 years certification of calibration, "buffered" probe
	in the product simulated solution, min/max memory. F/C switchable,
	field installable, and visual & audible temp alarm
	Standard Digital Temperature Display and Alarm Module with
	battery back-up to support temperature monitoring and alarms
	Sample simulated digital temperature display; .1° Resolution
	Adjustable temperature display in Celsius or Fahrenheit
	Audible and visual high/low temperature alarms
	2 probes; 1 in air and 1 in sample bottle
	Sensor error alarm
	Power failure alarm
	Resettable 1-point Min/Max History Remote alarm contacts
	Alarm validation (audible, visual, and remote alarm)
	Alarm mute with ring back
	5 total shelves, 3 adjustable plus 2 additional
	Vaccine storage power cord warning label
	Forced draft air circulation
	Keyed door lock
	Magnetic door gasket for positive seal
	Heavy duty door handle
	Swivel Casters (3 ½")
	LED interior lights are safety shielded and switch controlled
	Probe access port (3/8")
	Powder coated steel white interior and exterior
	Hydrocarbon, natural refrigerant (R290)
	Pharmacy refrigerator/freezer toolkit and temperature logs
	Two year parts and labor warranty, plus an additional three year
	compressor parts warranty
	Exterior dimensions: 25" W x 29 3/4" D x 65 3/4" H with casters; 62 1/4"
	H without casters
	Cycle defrost
	115V, 60 Hz, 3 Amps, 1/5 HP
	ETL listed
	Shipping weight: 258 lbs