

Serological Pipettes

Serological pipettes can be used to liquid transfer purpose of both tissue culture and bacterial culture experiments.

Features

1. Made of high clarity polystyrene.
2. Sterilized by E-beam
3. Non-pyrogenic, non-cytotoxic and non-hemolytic
4. DNase/RNase free
5. Polyolefin fiber filter to reduce contamination
6. Edged color-coded band for identification
7. Negative graduations for extra volume

Product Range



Individually Wrapped

Cat.No.	Description	Qty/Pk	Qty/Cs
DS324001	1ml Serological Pipette	500	3000
DS325001	2ml Serological Pipette	400	2400
326001	5ml Serological Pipette	200	800
327001	10ml Serological Pipette	200	800
328001	25ml Serological Pipette	200	800
DS329001	50ml Serological Pipette	100	600

Bulk, 25 Per Bulk Bag



Cat.No.	Description	Qty/Pk	Qty/Cs
DS326002	5ml Serological Pipette	500	1500
DS327002	10ml Serological Pipette	500	1500
DS328002	25ml Serological Pipette	200	800
DS329002	50ml Serological Pipette	100	600

Technical Drawing



