


THERMOTE[®] X NIST Traceable Wireless Temperature Sensor

Stay compliant with codes and regulations and reduce the chance of error or spoilage with the Thermote-X: A NIST traceable wireless ultra-low freezer sensor. The temperature range measures -328°F to 392°F (-200°C to 200°C) using a 78" probe. Get text alerts sent straight to your phone when your temperature goes out of the parameters you set. It has the option to connect via Wi-Fi or long-range LoRaWAN with a gateway (Sold separately).

The Thermote-X is durable. Rated IP67 our remote thermometers are water-proof and dust-proof. The Thermote-X is powered by 3 AA batteries (included). Unlike our competitors, our temperature monitor has low battery alerts. Battery lifespan depends on how often you receive data. Each sensor has a reset magnet that puts it back into configuration mode.

Never lose your records. Data logger automatically uploads data to the Zynect app or web portal when back on Wi-Fi. Create a unique code to share temperature sensor data and monitoring with friends, family, and coworkers. Zynect App available for iOS and Android. Compatible with IFTTT. No subscription. All sensors are made in Ithaca, NY.



Cat. #	Description	List
ZNX2500	Thermote [®] -X	\$319.00
ZNX2520	LoRaWAN Base Station	\$125.00
ZNXH1	Thermote [®] -X Mounting Holder	\$ 17.50