

TECHNICAL INFORMATION

dotBio

CRYOGENIC VIALS

External Thread Design



Internal Thread Design



Color Coded Cap Inserts




- Temps as low as -196°C
- Self-standing star foot
- Leakproof silicone O-ring
- Resealable bags of 50
- Gamma sterilized
- Free locking rack
- RNase, DNase, pyrogen, & endotoxin free
- Autoclavable

SPECIFICATIONS

Material (tube/cap):	LN ₂ grade polypropylene
Material (o-ring/washer):	Silicone
Temp Range:	-196°C to +121°C*
Vial Diameter:	12.5mm

*LN₂ gas phase only; not for liquid phase

Star foot base is compatible with locking rack



CRYOGENIC VIALS

Cat. No	Description	Cap Thread	Height (with cap)	Vol.	Qty	List Price
V3801	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	40mm	1mL	500/cs	\$131.00
V3802	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	47.3mm	2mL	500/cs	\$138.00
V3803	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	66.6mm	3mL	500/cs	\$175.50
V3805	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	External	90.7mm	5mL	500/cs	\$198.00
V4801	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	40.2mm	1mL	500/cs	\$131.00
V4802	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	47.5mm	2mL	500/cs	\$138.00
V4803	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	60.5mm	3mL	500/cs	\$175.50
V4805	Cryogenic Vials, Self-Standing, Sterile, Printed Graduations and Marking Area, Silicone O-Ring, 10 Bags of 50 Tubes	Internal	91mm	5mL	500/cs	\$198.00
V5809-A*	Color Cap inserts for Cryogenic Vials, assorted colors (white, blue, green, red, yellow), 500/pk	—	—	—	500/pk	\$21.00
V5920	Locking Rack for 50 "star foot" vials, 210 x 100 x 25mm	—	—	—	4/pk	\$49.00

ACCESSORIES

Cat. No	Description	Qty	List Price
R-2100	Freezer Box, hinged lid, polypropylene, 100 x 2 ml, 130 x 30 x 53 mm	5/pk	\$59.88
R2781-A*	FlipTop™ Freezer Box, hinged lid, cardboard, 81 x 2mL, 133 x 133 x 53mm	5/pk	\$17.70
R2700-A*	FlipTop™ Freezer Box, hinged lid, cardboard, 100 x 2mL, 133 x 133 x 53mm	5/pk	\$18.42

* Assorted color pack, also available in individual colors

FREE
Locking Rack
in every case

dot scientific inc

www.dotscientific.com

(800) 878-1785

orders@dotscientific.com