



Levo Plus

Features

- Intuitive and convenient speed with 6 levels
- Lightweight pipette filler compatible with glass or plastic serological pipettes (0.1 – 100ml)
- Ergonomic design for fatigue-free and effective streamlined pipetting
- LCD screen displays battery status and pipetting speed
- Quiet and powerful motor and pump fills 25 mL pipette in less than 5 seconds
- Long-lasting lithium-ion battery

Levo Plus is going green

- New package size is 65% smaller to help reduce the carbon footprint of shipping.
- Optional Accessories can be purchased separately: Wall-mounted Holder; Standing Holder and Hydrophobic Filter



Special offer \$99

Part number: 7033110300
Promo ends 12.31.2022



dot scientific inc.

4165 Lippincott, Burton, MI 48519
800-878-1785 | orders@dotscientific.com