

## WORKSPACE SERIES

Own your workspace and its efficiency. From a busy lab to a complicated control panel, your work requires products that fit into your space, not the other way around. This workspace series makes room for two different high-performing solutions.

### M210-LAB Label Printer

The M210-LAB label printer takes on your toughest labeling applications in the lab. Whether you're labeling beakers for short experiments or dropping vials in a cryo freezer, the label materials and built-in label creation guides are a perfect partnership for your laboratory. And with smart automatic formatting, you're set to go in seconds.



#### Top applications



##### Laboratory

Vial, tube, slide labeling, sample identification



##### General identification

Directional marking, equipment labeling

#### Unique features

- Automatic label formatting, including five pre-loaded application design guides, for rapid and accurate lab label creation
- Ability to store 12 label files up to 2MB and an on-board library of 133 symbols
- Hot keys and menu options for fast label creation
- Ability to create UPC, DataMatrix and QR Codes
- 90+ compatible easy-to-load cartridges (continuous and pre-sized), including materials specifically designed to withstand lab environments

#### Label output

**Speed:** Prints up to 100 labels / day at 0.4 inches (10.2 mm) per second

**Quality:** 203 dpi print resolution

#### Optional accessories

- Multi-functional tool that includes a heavy-duty magnet, retractable stand and flashlight
- Attachable magnet for hands-free label creation
- Rechargeable Li-Ion battery pack
- Hard and soft carrying cases to protect your label printer and keep accessories organized

#### Available kits

Printer with accessory kit



**scientific inc.**

4165 Lippincott, Burton, MI 48519  
800-878-1785 | [orders@dotscientific.com](mailto:orders@dotscientific.com)