Product Name:	PH-ABT-HC-10PG
Description:	Premier Pharmacy/Vaccine Compact Refrigerator 10.5 Cu. Ft. Capacity
	One glass swing door, right hinged, not reversible
	Microprocessor temperature controller
	Digital temperature display
	Audible and visual high and low temperature alarms
	2 probes; 1 in air and 1 in sample bottle
	Remote alarm contacts
	Adjustable operating temperature range: 36°F – 46°F (2°C -8°C)
	Temperature monitor device (TMD) complies with the current CDC
	guidelines, with 3 years certification of calibration, "buffered" probe
	in the product simulated solution, min/max memory. F/C switchable,
	field installable, and visual & audible temp alarm
	Six adjustable shelves with guard rail on back to prevent product
	from being pushed against back wall
	Bottom footed shelf to allow air circulation
	Back wall evaporator cover
	Vaccine storage power cord warning label
	Forced draft circulation
	Keyed door lock
	Magnetic door gasket for positive seal Leveling legs
	LED interior light
	Probe access port (3/8")
	Powder coat exterior
	High density urethane foam cabinet insulation
	HFC-free refrigerant (R600a), foam insulation and packaging
	Pharmacy refrigerator/freezer toolkit and temperature logs
	Two year parts and labor warranty, plus an additional three year
	compressor parts warranty
	Pyxis [®] , Omnicell [®] and AcuDose RX [®] compatible
	Exterior dimensions: 23 ³ / ₄ "W x 26 ¹ / ₂ "D x 59"H
	Note: This unit must have 4" clearance on sides and back for
	adequate ventilation
	UL/C-UL listed
	ENERGY STAR CERTIFIED
	Cycle defrost
	115V, 60 Hz, 1/3 HP, 2.3 Amps
	Energy Consumption 2.42 kWh/day
	Shipping weight: 209 lbs.