Product Name:	ABT-FRP-04
Description:	Standard Undercounter Flammable Storage Refrigerator Freestanding 5 Cu. Ft. capacity One solid swing door, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Adjustable operating temperature range: 1°-10°C Remote alarm contacts Three adjustable shelves with guard rail on back to prevent product from being pushed against the back wall Keyed door lock Magnetic door gasket for positive seal Integrated door handle One probe access port (3/8") White cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty Note: This unit must have 4" clearance on sides and back for adequate ventilation Exterior dimensions: 23 5/8" W x 24 ¼ " D x 32 3/16" H Cycle defrost 115V, 60 Hz, 1.2 Amps, 1/9 HP UL/C-UL listed Units meet NFPA guidelines 45 and 70 Shipping weight: 118 lbs.