

Product Name:	ABT-FFP-04
Description:	<p>Standard Undercounter Flammable Storage Freezer Freestanding 4 Cu. Ft. capacity One solid swing door, right hinged, not reversible Microprocessor temperature controller Digital temperature display Audible and visual high and low temperature alarms 2 probes; 1 in air and 1 in sample bottle Remote alarm contacts Adjustable operating temperature range: -15°-25°C 2 fixed shelves with three removable baskets Keyed door lock Magnetic door gasket for positive seal Integrated door handle One probe access port (3/8") White cabinet and door High density urethane foam cabinet and door insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> Exterior dimensions: 23 5/8"W x 24 1/4"D x 32 3/16" H Manual defrost 115V, 60 Hz, 1.4 Amps, 1/8 HP UL/C-UL listed Units meet NFPA guidelines 45 and 70 Shipping weight: 118 lbs.</p>