

Product Name:	ABT-HC-10PG
Description:	<p>Premier Glass Door Compact Laboratory Refrigerator 10.5 Cu. Ft. Capacity One glass swing door, right hinged, not reversible Microprocessor temperature controller Digital temperature display with sample simulated probe Audible and visual high and low temperature alarms Remote alarm contacts Adjustable operating temperature range: 1°C to 10°C Six adjustable shelves with guard rail on back to prevent product from being pushed against back wall Bottom footed shelf to allow air circulation Forced draft air circulation Back wall evaporator cover Keyed door lock Magnetic door gasket for positive seal Leveling legs LED interior light Probe access port (3/8") Powder coat exterior High density urethane foam cabinet insulation HFC-free refrigerant (R600a), foam insulation and packaging Two year parts and labor warranty, plus an additional three year compressor parts warranty Exterior dimensions: 23 ¾"W x 26 ½"D x 59"H <i>Note: This unit must have 4" clearance on sides and back for adequate ventilation</i> UL/C-UL listed ENERGY STAR CERTIFIED Cycle defrost 115V, 60 Hz, 1/3 HP, 2.3 Amps Energy Consumption 2.42 kWh/day Shipping weight: 209 lbs.</p>