



Key Features

- Microprocessor temperature controller
- Digital temperature display
- Audible and visual high and low temperature alarms
- Remote alarm contacts
- 2 probes; 1 in air and 1 in sample simulated bottle
- Adjustable operating temperature range:
 - Refrigerator: 36°F – 46°F (2°C -8°C)
 - Freezer: 4°F – -13°F (-15°C - -25°C)
- Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable
- Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage
- Keyed door lock
- Magnetic door gasket for positive seal
- Leveling legs
- Probe access port (3/8")
- High density urethane foam cabinet and door insulation
- Hydrocarbon natural refrigerant (R600a)
- Pharmacy refrigerator/freezer toolkit and temperature logs
- Pyxis®, Omnicell® and AcuDose RX® compatible
- Quality certified with 100% unit temperature testing
- Defrost:
 - Refrigerator: Cycle defrost
 - Freezer: Auto & Manual defrost available
- ADA compliant units available
- Two year parts and labor warranty
- UL/C-UL listed

BUILT-IN UNDER COUNTER REFRIGERATORS AND FREEZERS

Front Ventilation design for built-in applications between pharmacy cabinetry




The Pharmacy Premier Built-In Under Counter Refrigerators & Freezers take a large step forward in performance from past solid door general purpose units; offering a digital temperature display and alarm module coupled with a microprocessor temperature controller.

An interior fan driven, forced draft air circulation system provides excellent temperature uniformity and rapid temperature recovery after door openings.

Rear wall evaporator cover.



Pharmacy Premier Built-In Under Counter Refrigerators & Freezers



Now Featuring!
Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration



Premier Refrigerator

Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
PH-ABT-HC-UCBI-0404	4.6	Cycle	1 Solid	2 Adj/1 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.74
PH-ABT-HC-UCBI-0404-LH	4.6	Cycle	1 Solid	2 Adj/1 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.74
PH-ABT-HC-UCBI-0404G	4.6	Cycle	1 Glass	2 Adj/1 Fixed	23 3/4	25 1/2	33 3/8	R600a	1/8	1.74
PH-ABT-HC-UCBI-0404G-LH	4.6	Cycle	1 Glass	2 Adj/1 Fixed	23 3/4	25 1/2	33 3/8	R600a	1/8	1.74
PH-ABT-HC-UCBI-0404G-ADA	4.6	Cycle	1 Glass	2 Adj/1 Fixed	23 3/4	25 1/2	31 15/16	R600a	1/8	1.74
PH-ABT-HC-UCBI-0404G-ADA-LH	4.6	Cycle	1 Glass	2 Adj/1 Fixed	23 3/4	25 1/2	31 15/16	R600a	1/8	1.74
PH-ABT-HC-UCBI-0404SS	4.5	Cycle	1 Solid	2 Adj/1 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.74
PH-ABT-HC-UCBI-0404SS-LH	4.5	Cycle	1 Solid	2 Adj/1 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.74

Premier Freezer


Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
PH-ABT-HC-UCBI-0420	4.2	Manual	1 Solid	2 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.5
PH-ABT-HC-UCBI-0420-LH	4.2	Manual	1 Solid	2 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.5
PH-ABT-HC-UCBI-0420-ADA	4.2	Manual	1 Solid	2 Fixed	23 3/4	24 1/2	32	R600a	1/8	1.5
PH-ABT-HC-UCBI-0420-ADA-LH	4.2	Manual	1 Solid	2 Fixed	23 3/4	24 1/2	32	R600a	1/8	1.5
PH-ABT-HC-UCBI-0420SS	4.2	Manual	1 Solid	2 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.5
PH-ABT-HC-UCBI-0420SS-LH	4.2	Manual	1 Solid	2 Fixed	23 3/4	24 1/2	33 3/8	R600a	1/8	1.5
PH-ABT-HC-UCBI-0420A	4.2	Auto	1 Solid	3 Adj	23 3/4	24	33 1/8	R600a	1/8	1.8
PH-ABT-HC-UCBI-0420A-LH	4.2	Auto	1 Solid	3 Adj	23 3/4	24	33 1/8	R600a	1/8	1.8

Pharmacy Premier Freestanding Under Counter Refrigerators & Freezers



Standard features include a microprocessor temperature controller and display for superior temperature control, verification, and recovery.

Now Featuring!
Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration



Digital temperature controller allows for very precise temperature maintenance

Certain models available in either glass or solid door

PH-ABT-HC-UCFS-0204



Premier Refrigerator

Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
PH-ABT-HC-UCFS-0104	1	Cycle	1 Solid	2 Adj	17 1/4	19 1/4	21 1/4	R600a	1/8	0.9
PH-ABT-HC-UCFS-0104-LH	1	Cycle	1 Solid	2 Adj	17 1/4	19 1/4	21 1/4	R600a	1/8	0.9
PH-ABT-HC-UCFS-0104G	1	Cycle	1 Glass	2 Adj	17 1/4	19 1/4	21 1/4	R600a	1/8	0.9
PH-ABT-HC-UCFS-0104G-LH	1	Cycle	1 Glass	2 Adj	17 1/4	19 1/4	21 1/4	R600a	1/8	0.9
PH-ABT-HC-UCFS-0204	2.5	Cycle	1 Solid	4 Adj	17 1/4	19 1/4	28 7/8	R600a	1/8	0.9
PH-ABT-HC-UCFS-0204G	2.5	Cycle	1 Glass	4 Adj	17 1/4	19 1/4	28 7/8	R600a	1/8	0.9
PH-ABT-HC-UCFS-0504G	5.2	Cycle	1 Glass	3 Adj	23 3/4	24 1/2	32 1/8	R600a	1/6	1.3

Premier Freezer

Model	Cu. Ft.	Defrost	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
PH-ABT-HC-UCFS-0120	1.7	Manual	1 Solid	2 Adj	23 3/4	24	20 3/8	R600a	1/8	1
PH-ABT-HC-UCFS-0120-LH	1.7	Manual	1 Solid	2 Adj	23 3/4	24	20 3/8	R600a	1/8	1
PH-ABT-HC-UCFS-0120A	1.3	Auto	1 Solid	2 Adj	23 1/2	23 1/2	20 1/2	R600a	1/7	1.6
PH-ABT-HC-UCFS-0120A-LH	1.3	Auto	1 Solid	2 Adj	23 1/2	23 1/2	20 1/2	R600a	1/7	1.6

Key Features

Freestanding UNDER COUNTER & COUNTER TOP REFRIGERATORS AND FREEZERS

The ideal choice for pharmacies with under counter and counter top applications with adequate space for ventilation



PH-ABT-HC-UCFS-0504G

An interior fan driven, forced draft air circulation system provides excellent temperature uniformity and rapid temperature recovery after door openings.



PH-ABT-HC-UCFS-0104



PH-ABT-HC-UCFS-0104G

- Microprocessor temperature controller)
- Digital temperature display
- Audible and visual high and low temperature alarms, with battery back up (batteries included)
- Remote alarm contacts
- Adjustable operating temperature range:
 - Refrigerator: 36°F – 46°F (2°C -8°C)
 - Freezer: 4°F – -13°F (-15°C - -25°C)
- Refrigerator/Freezer freestanding thermometer with 3 year certificate of calibration, with min/max memory, °F/°C switchable
- Cord label stating Warning! Do not unplug the appliance or break circuit. Expensive vaccine in storage
- Keyed door lock
- Back wall evaporator cover
- Forced draft circulation (refrigerators only)
- Magnetic door gasket for positive seal
- Leveling legs
- Probe access port (3/8")
- Powder coat exterior
- High density urethane foam cabinet and door insulation
- Hydrocarbon natural refrigerant (R600a)
- Pharmacy refrigerator/freezer toolkit and temperature logs
- Pyxis®, Omnicell® and AcuDose RX® compatible
- Quality certified with 100% unit temperature testing
- Note: This unit must have 4" clearance on sides and back for adequate ventilation
- Defrost:
 - Refrigerator: Cycle defrost
 - Freezer: Auto & Manual defrost available
- ADA compliant units available
- Two year parts and labor warranty
- UL/C-UL listed

Requires at least 4" of clearance on each side for adequate ventilation