



the belly button

from



specifications:

● Base Footprint:	8"x8"
● Standard Platform Size	7"x7"
● Overall Height	5.25"
● Product Weight	5lbs
● Shipping Weight	7lbs
● Maximum Weight on Platform	5lbs
● Electrical (Universal Input)	94-240V
● Motor (Brushless DC)	12VDC
● Operating Temperature Range	0°C - 60°C
● Maximum Operating Humidity (Non-Condensing)	97%
● Operates in CO ₂ Incubator	Yes
● Operational Speed	0 - 100RPM
● Platform Pitch Control	0° - 8°



features:

- Small footprint takes up little lab space and easily fits inside incubators.
- Flat 115V power cord fits through incubator doors without compromising environment.
- Unique undulating motion that can be adjusted with the pitch control collar from 0 to 8 degrees.
- Variable speed with smooth operation from 6 to 100RPM, and no cogging at slower speeds.
- Comes with white, non-skid mat that can be easily removed and cleaned.
- The wide operational temperature range of 0° to 60°C allows the unit to be used in high-humidity incubators and cold rooms.
- High-Quality 12VDC brushless DC motor provides cool operation and years of dependable use.
- Platform is made of chemically resistant 3/8" cast acrylic with no rusting or discoloration, like metal-top units.
- Finger Mat is available for shaking 10ml to 50ml tubes.

applications:

- Cell Culture Growth/Mixing
- Immunoprecipitation
- Northern and Southern Hybridizations
- Phage Elution and Amplification of cDNA Libraries
- Western Blotting and Membrane Rinsing
- Staining and De-Staining Gels
- Miscellaneous Mixing

ordering information:

CATALOG#	DESCRIPTION
BBUAAUV1S	THE BELLY BUTTON - UV, 50/60Hz (CE)
BBULS0001	THE BELLY BUTTON, w/TUBE MAT - UV, 50/60Hz (CE)



4165 LIPPINCOTT ■ BURTON, MI 48519

TOLL FREE: 1-800-878-1785 ■ FAX: 1-810-744-2865

E-MAIL: sales@dotscientific.com ■ WEB: www.dotscientific.com