



S2031-13 with K566



Specifications

Temp Range	Ambient +5° to 80°C
Temp Uniformity	±0.25°C
Temp Accuracy	±0.1°C
Shaker speed	20 - 300 rpm
Timer	0-99 min or continuous
Shaker orbit	19 mm
Exterior dims	42.5 x 55 x 58 cm/ 16.75 x 21.75 x 23 in
Chamber dims	34.3 x 37.5 x 38.1 cm/ 13.5 x 14.25 x 15 in
Chamber volume	1.7 cu ft/49 liters
Weight	83 lb/37.7 kg



- Integrated digital orbital shaker
- Exclusive SmartCheck temperature control system
- Broad temperature ranges to meet a wide variety of applications
- Stacking design saves space

Ordering Information

Cat#	Description
I5211-DS	211DS Shaking Incubator, 120V
I5211-SA	Stacking adapter, required when stacking incubators
I5230	Flask clamp platform, pre drilled (see clamps below)
I5231	Flat platform with non-slip rubber mat (30 x 30cm)
I5231-D	Double flat platform with non-slip rubber mat
S2031-13	Universal spring loaded clamp platform for flasks, bottles, etc.
S2040-85	Flask clamp, 125mL (max 16)
S2040-99	Flask clamp, 250mL (max 9)
S2040-09	Flask clamp, 500mL (max 5)
S2040-01	Flask clamp, 1000mL (max 4)

* To order 230V units add -230V to the end of the catalog number. 230V unit includes -EU and -UK cords
 **Other size flask clamps are available. Call for details.