

M48725Q Vari-Mix™ and Speci-Mix™ Test Tube Rockers

Obtain uniform specimen suspension with the smooth rocking action of Thermo Scientific™ Vari-Mix™ and Speci-Mix™ Test Tube Rockers. Compact and precision-controlled platform mixers, all models hold 10 to 30mm diameter tubes and are ideal for blood collection tubes.

Product Size	Each
Model	Vari-Mix Test Tube Rocker
Tube Capacity	16 (10-20mm); 7 (22-30mm)
Speed Range	2 to 20rpm
Temperature Range (Metric)	17° to 27°C (Operating)
Voltage	120V
Rocking Angle	Adjustable rocking from 1° to 48° angle
Humidity Range	20 to 80%
Height (English)	4.3 in.
Length (English)	5.5 in.
Width (English)	15 in.
Height (Metric)	10.7cm
Length (Metric)	13.9cm
Width (Metric)	38.1cm
Shipping Weight (English)	5.8 lb.
Shipping Weight (Metric)	2.6kg
Hertz	50Hz
Amperage	0.08A
Dimensions (L x W x H)	5.5 x 15 x 4.3 in. (13.9 x 38.1 x 10.7cm)
Electrical Requirements	120V 50Hz
Speed	2 to 20rpm
Type	Rockers