


NEW! **VORTEX-GENIE® PULSE** Pulsing Vortex Mixer



Vortex-Genie Pulse Accessories

More Than Just Vortexing!

Another revolutionary product from the makers of the Vortex-Genie!

- ⦿ True pulsing action provides more collisions and random motion.
- ⦿ For applications requiring more aggressive action than standard vortexing.
- ⦿ Variable "RUN/STOP" time provides quick pulsing or interval mixing.
- ⦿ Fully programmable "ON" and "REST" function reduces heat buildup within the vessel, while keeping materials in suspension. Ideal for heat sensitive samples.
- ⦿ Digital control and display of speed, time, and "RUN/STOP" time, for reproducibility and accuracy, with true speed display.
- ⦿ Use with all Vortex-Genie accessories.
- ⦿ May also be used as a standard "continuous" Vortex Mixer.
- ⦿ Ideal for applications such as re-suspensions of lyophilized materials, creating solutions, enzymatic reactions and much more.

Ordering Information

Catalog No.	Description
SI-P236	Vortex-Genie Pulse, 120V
SI-P246	Vortex-Genie Pulse, 230V, No Plug
SI-P256	Vortex-Genie Pulse, 230V, Euro Plug
SI-P266	Vortex-Genie Pulse, 230V, British Plug
SI-P276	Vortex-Genie Pulse, 230V, Swiss Plug
SI-P296	Vortex-Genie Pulse, 240V, Australian Plug
SI-P286	Vortex-Genie Pulse, 100V

Specifications

Supplied	3-inch Platform, Pop-off Cup
Speed	500 - 3000 RPM (2850 for 50Hz models)
Timer - Hands Free	1-99 minutes or continuous
Pulse - ON/OFF	0.1 seconds - 99 minutes
Weight	4 kgs (8.8 lbs)
Dimensions (D x W x H)	165 x 122 x 165mm (6.5 x 4.8 x 6.5in)



Home of the Multi-Task Genie™ Family
 70 Orville Drive, Bohemia, NY 11716 U.S.A.
 Phone: (631) 567-4700 • Fax: (631) 567-5896
 Toll Free: 888-850-6208
 Customer Service: custsvc@scientificindustries.com
scientificindustries.com

Multi-Task Genie™ Family Accessories

A Vast Range Of Accessories For The Ultimate Multi-Task Versatility In Any Application!



Accessories for all Vortex-Genie® 2 mixers and Vortex-Genie® Pulse

Fig. 1
Model H301
Multiple Sample
Starter Set




Fig. 4
Model H302
Large Sample Set




Fig.	Cat. No.	Description	Fig.	Cat. No.	Description
1	504-0233-00	Model H301 Multiple Sample Starter Set	--	580-2013-00	Rubber Cover for 3-inch Platform (not shown)
--	146-6005-00	6-inch Platform (not shown separately)	13	SI-0564	TurboMix™ Attachment
2	504-0234-00	Microtube Foam Insert (2)**	14	SI-0565	Universal Microtube Holder*
3	504-0235-00	Microplate Foam Insert (2)**	15	SI-0511	Large Ampule/Tube Attachment
4	503-0282-00	Model H302 Large Sample Set	16	SI-0570	Small Ampule/Tube Attachment
5	504-0039-00	Recessed Platform	--	SI-0510	Two-Tier Microplate Foam Insert (not shown)
6	503-0278-00	9-16mm Tube Foam Inserts (2)***	17	SI-H506	Horizontal 50ml Tube Holder
7	503-0279-00	14-29mm Tube Foam Inserts (2)***	18	SI-H512	Horizontal 15ml Tube Holder
8	503-0280-00	29-37mm Tube Foam Inserts (2)***	19	SI-H524	Horizontal Microtube Holder
9	503-0277-00	Blank Foam Inserts (2)***	20	SI-V506	Vertical 50ml Tube Holder
--	503-0283-00	Assorted Foam Inserts (4)***	21	SI-V512	Vertical 15ml Tube Holder
10	568-0001-00	Elastic Bands (2)	22	SI-V524	Vertical Microtube Holder
11	146-3011-00	Pop-off Cup	23	SI-V525	Vertical 50ml, 15ml, Microtube Holder
12	0K-0500-902	3-inch Platform with Rubber Cover	--	0K-0236-408	Bottom Closure with Screws Kit

* May be used with TurboMix attachment or Disruptor Genie®

** 6-inch Platform Required

*** Recessed Platform Required

NOTE: Vortex-Genie 1 may only be used with Pop-off Cup and 3-inch Platform

U.S. Pat. 5,707,861



70 Orville Drive, Bohemia, NY 11716 U.S.A.
Phone: (631) 567-4700 • Fax: (631) 567-5896
Toll Free: 888-850-6208
Customer Service: custsvc@scientificindustries.com
scientificindustries.com