


Standard Manual Defrost Freezers

Standard Manual Defrost Freezers

The Standard Manual Defrost Freezers are intended for general purpose applications where precise temperature control and maintenance are not essential and there are no requirements for temperature alarms. Manual defrost units have been designed to operate without product damaging auto defrost cycles. **Periodic manual defrost is required on these units.**



Model	Cu. Ft.	Defrost	Temp	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
ABT-MFS-14	14	Manual	-20°C	1 Solid	3 Fixed	29 3/4	32 3/8	62 1/2	R134a	1/3	5
ABT-MFS-17	17	Manual	-20°C	1 Solid	3 Fixed	34	30 7/16	67 3/4	R134a	1/3	5
ABT-MFS-20	20	Manual	-20°C	1 Solid	4 Fixed	34	32 9/16	73 1/4	R134a	1/3	5
ABT-MFS-30	30	Manual	-20°C	1 Solid	3 Fixed/ 2 Adj.	35 3/4	34 1/4	72 1/2	R404a	1/3	6
Chest Freezers											
Model	Cu. Ft.	Defrost	Temp	Doors	Shelves	W"	D"	H"	RL	H.P.	Amps
ABT-MFS-16-C	16	Manual	-20°C	1 Hinge	1 Basket	56	29 1/2	35	R134a	1/4	5
ABT-MFS-22-C	22	Manual	-20°C	1 Hinge	4 Baskets	73 1/4	29 1/2	35	R134a	1/4	5

- One solid swing door, right hinged, not reversible (Upright Freezers)
- Mechanical thermostat
- Adjustable operating temperature range: -15° -25° C
- Keyed door lock
- Magnetic door gasket for positive seal
- Heavy duty door handle
- Leveling legs *Not included on Chest Freezer models*
- Probe access port (3/8")
- White textured cabinet and door
- High density urethane foam cabinet and door insulation
- CFC-free refrigerant (R134a), foam insulation and packaging
- One year parts and labor warranty; plus additional four year compressor parts warranty
- Quality certified with 100% unit temperature testing
- Manual defrost
- UL/C-UL listed