









Thermo Scientific TSE and TSD Series

Thermo Scientific TSE and TSD Series Upright Freezer Racks



2" Box Racks	Model No. (Description)	Dimensions H x W x D in. (cm)	Freezer Vol. cu. ft. (liters)	Boxes/Rack	Racks/Shelf	Racks/Freezer	Boxes/Freezer
	920094 (Sliding Drawer Rack for 2" Boxes)	11.9 x 5.5 x 16.6 (30.2 x 14 x 42.2)	13 (368.1)	15	4	16	240
	920096 (Sliding Drawer Rack for 2" Boxes)	11.9 x 5.5 x 22 (30.2 x 14 x 55.9)	17.3 (490)	20	4	16	320
			23 (651.3)	20	5	20	400
	920090 (Sliding Drawer Rack for 2" Boxes)	11.9 x 5.5 x 26.3 (30.2 x 14 x 66.5)	28.1 (793)	25	6	24	600
	920015 (Adjustable Side Access Rack)	10.9 x 5.5 x 16.5 (27.7 x 14 x 41.9)	13 (368.1)	15	4	16	240
	920022 (Side Access Rack for 2" Boxes)	10.9 x 5.5 x 24.6 (27.7 x 14 x 62.5)	17.3 (490)	22	4	16	352
			23 (651.3)	22	5	20	440
	920018 (Adjustable Side Access Rack)	10.9 x 5.5 x 22 (27.7 x 14 x 55.9)	17.3 (490)	20	4	16	320
			23 (651.3)	20	5	20	400

Microplate Racks	Model No. (Description)	Dimensions H x W x D in. (cm)	Freezer Vol. cu. ft. (liters)	Plates/Rack	Plates/Shelf	Plates/Freezer	Plates/Freezer
	398335 (Side Access Deepwell Rack)	9.5 x 5.5 x 22 (24 x 14 x 60)	17.3 (490)	30	4	16	480
			23 (651.3)	30	5	20	600
	820036 (Side Access Deepwell Rack)	12 x 5.5 x 22 (30.5 x 14 x 60)	17.3 (490)	36	4	16	576
			23 (651.3)	36	5	20	720
	820066 (Side Access Deepwell Rack)	9.5 x 5.5 x 22 (24 x 14 x 60)	17.3 (490)	66	4	16	1,056
			23 (651.3)	66	5	20	1,320

3" Box Racks	Model No. (Description)	Dimensions H x W x D in. (cm)	Freezer Vol. cu. ft. (liters)	Boxes/Rack	Racks/Shelf	Racks/Freezer	Boxes/Freezer
	920095 (Sliding Drawer Rack for 3" Boxes)	11.9 x 5.5 x 16.5 (30.2 x 14 x 41.9)	13 (368.1)	9	4	16	144
	920097 (Sliding Drawer Rack for 3" Boxes)	11.9 x 5.5 x 22 (30.2 x 14 x 55.9)	17.3 (490)	12	4	16	192
			23 (651.3)	12	5	20	240
	920091 (Sliding Drawer Rack for 3" Boxes)	11.9 x 5.5 x 26.3 (30.2 x 14 x 66.5)	28.1 (793)	15	6	24	360
	920009 (Side Access Rack for 3" Boxes)	10.9 x 5.5 x 16.5 (27.7 x 14 x 41.9)	13 (368.1)	9	4	16	144
	920017 (Side Access Rack for 3" Boxes)	10.9 x 5.5 x 22 (27.7 x 14 x 55.9)	17.3 (490)	12	4	16	192
			23 (651.3)	12	5	20	240

Sliding Drawer Inserts (Pack of 5)

Model No.	For	Freezer Volume	Slots	Model No.	For	Freezer Volume	Slots
1950666P5	microplates	13 cu. ft. (368L)	4	1950668P5	boxes	13 cu. ft. (368L)	3
1950669P5	microplates	17/23 cu. ft. (490/651L)	6	1950670P5	boxes	17/23 cu. ft. (490/651L)	4
1950671P5	microplates	28 cu. ft. (793L)	7	1950672P5	boxes	28 cu. ft. (793L)	5

Specialty Racks	Model No. (Description)	Dimensions H x W x D in. (cm)	Freezer Vol. cu. ft. (liters)	Racks/Shelf	Racks/Freezer
	830012 (Multi-Drawer Rack)	11.5 x 11 x 11 (29.2 x 28 x 28)	13 (368.1)	1	4
			17.3 (490)	2	8
			23 (651.3)	4	16
	830030 (8 Compartment Container)	5.3 x 7 x 18 (13.5 x 18 x 45.7)	13 (368.1)	6	24
	830032 (10 Compartment Container)	5.3 x 7 x 22.6 (13.5 x 18 x 56)	17.3 (490)	6	24
			23 (651.3)	6	32